



WHEATON
CRIMPenstein™
Electric Crimper/Decapper





WHEATON CRIMPenstein™ Electric Crimper/Decapper

**Advanced Ergonomic Design
Engineered to Save Money
and Make Your Crimping and
Decapping Jobs Go Faster
...and Easier**

Lightweight and easy to operate – Enables handheld semi-automatic operation for extended periods with less effort

Quick-change jaw-sets available in 8, 11, 13, and 20 mm sizes

– Provides cost savings by eliminating the need to purchase separate crimpers and decappers

Electric powered for quick setup and high volume operation

– Eliminates the need for regulated compressed air or batteries which discharge over time

Adjustable Jaw Intensity with Memory – Allows user to regulate jaw action as required to achieve the perfect crimp

Electrify your crimping and decapping operations with the new innovative patent pending Wheaton Crimpenstein. Engineered to improve productivity and efficiency in a variety of laboratories performing chromatography, QC testing, lyophilization and vaccine development, or any other applications where vials are capped with aluminum seals. This AC powered crimper/decapper provides fast and consistent results with a simple push of a button. Jaw action intensity can be adjusted as needed. The ability to change jaw-sets offers an unmatched level of flexibility and economy for applications where various sizes of aluminum seals are required. A sturdy stand is available which will hold the Crimpenstein in a fixed vertical position for operations where vials need to be handled one at a time.

- Interchangeable jaw-sets save time and money for multiple application laboratories

- Ideal for high volume handheld operation, the lightweight Crimpenstein is ergonomically designed with your comfort and safety in mind



- The crimper stand holds the Crimpenstein in a fixed vertical position enabling individual capping to be performed quickly and easily



Ordering Information

Crimpenstein™ Unit and Power Cord*

Description	Size (mm)	Cat. No.	Qty/Cs
With Crimper and Decapper Jaw-Set			
Crimper and Decapper	8	W225600-A	1
Crimper and Decapper	11	W225601-A	1
Crimper and Decapper	13	W225602-A	1
Crimper and Decapper	20	W225603-A	1
With One Crimper Jaw			
Crimper	8	W225605-A	1
Crimper	11	W225606-A	1
Crimper	13	W225607-A	1
Crimper	20	W225608-A	1
With One Decapper Jaw			
Decapper	8	W225610-A	1
Decapper	11	W225611-A	1
Decapper	13	W225612-A	1
Decapper	20	W225613-A	1

*Detachable power cord

Crimper and Decapper Jaws Sold Separately**

Description	Size (mm)	Cat. No.	Qty/Cs
Crimper and Decapper Jaw-Set			
Crimper and Decapper	8	W225615	1
Crimper and Decapper	11	W225616	1
Crimper and Decapper	13	W225617	1
Crimper and Decapper	20	W225618	1
Crimper Jaw Sold Separately			
Crimper Jaw	8	W225620	1
Crimper Jaw	11	W225621	1
Crimper Jaw	13	W225622	1
Crimper Jaw	20	W225623	1
Decapper Jaw Sold Separately			
Decapper Jaw	8	W225625	1
Decapper Jaw	11	W225626	1
Decapper Jaw	13	W225627	1
Decapper Jaw	20	W225628	1

**Not a complete unit



Each Crimpenstein comes with a convenient carrying/storage case



Crimpenstein Stand

Description	Cat. No.	Qty/Cs
Crimpenstein Stand Only	W225630	1

For a complete list of serum bottles, vials and closures see our on-line catalog...

www.wheatonsci.com



WHEATON

Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038

1.800.225.1437 ■ 1.856.825.1100 ■ 1.856.825.1368 (F)

www.wheatonsci.com